

FarrowWrap® Trim-To-Fit Knee High Order Form



Date://	Patient Name:	Therapist Name:
Ship to:		
R: L: Measure Posteriorly leg contour	B least ankle	R. L. B circumference
	Straight distance *	R: L: X straight distance

FarrowWrap® TTF Knee High

FarrowWrap® TTF Legpiece		Medium		XLarge		FarrowWrap® TTF Footpiece Mediu	Medium	
(L4) & (L5) Circumferences (L3) Circumference		20 – 60 cm 45 – 90 cm 20 – 55 cm 45 – 90 cm			Forefoot (X) circumference $20-30$ (X) Heel to forefoot length $14-21$			
(L2) Circumference(B) Ankle Circumference2B-D) Posterior Leg Length		20 – 45 cm 20 – 40 cm 27 – 35 cm		35 – 80 cm 35 – 70 cm 27 – 35 cm		LITE TTF Footpiece FWLT-T-FT3 Qty STRONG TTF Footpiece FWST-T-FT3 Qty		
LITE TTF Legpiece STRONG TTF Legpiece	FWLT-T-LT3	Qty Qty		Qty		Please note that you may substitute a FarrowHybrid Foot Compression for a footpiece. See FarrowHybrid order form for additional details.		

FarrowWrap® TTF garments fit a very wide range of sizes and can be trimmed on site for a perfect fit, and retrimmed if the patient's edema further reduces.

When measuring the circumferences from B to L5, the distance between each circumference will vary for the different fabrics. For LITE legpieces, measure the circumference at **point B** and then measure the circumference every **6.5 cm** until measurement **point D** is reached. For **STRONG** legpieces, measure the circumference at **point B** and then measure the circumference every **5 cm** until measurement **point D** is reached/

FarrowWrap® TTF (Trim-To-Fit) garments require trimming by a professional before being applied by the patient. When a TTF footpiece and legpiece are purchased together, a free pair of liners is included. Silver liners are supplied by default, but you may request Terry Cloth instead.