



Date: \_\_\_/\_\_\_/\_\_\_ Patient Name: \_\_\_\_\_ Therapist Name: \_\_\_\_\_

Ship to: \_\_\_\_\_

**Measure in centimeters**

Left  Right

**Length Measurement from axilla to wrist (c-g) = \_\_\_\_\_**

**LONG LENGTH IS > 43 CMS      REGULAR LENGTH IS < 43 CMS\***

The length measurement is taken from **Wrist** along the inside of the arm to the inseam, **Axilla**, of the arm; measure with patient standing or sitting with arm down by the side.

Measurement Points

Circumference **Axilla** is taken around the upper arm at the axilla. For ease in determining this point, place a book up into the arm pit area. Point 'g' on the arm will be even with the top of the book.

Circumference **Wrist** is the point of greatest compression and therefore the most important point. This point is found distally on the arm where the wrist meets the hand. One should find the smallest circumference at this point. **Be sure not to measure too tight.**

**Circumference Wrist** = \_\_\_\_\_ with Silicone\*  with no Silicone

**Circumference 'Mid-Lower'** = \_\_\_\_\_

**Circumference 'Mid-Upper'** = \_\_\_\_\_

**Circumference Axilla** = \_\_\_\_\_

Material	Class I 15-20 mmHg	Class II 20-30 mmHg	Class II 30-40 mmHg	Colors Available
Bella Strong	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/>
Bella Lite Arm Sleeve	<input type="checkbox"/>	<input type="checkbox"/>		Beige <input type="checkbox"/> Brown <input type="checkbox"/> Black <input type="checkbox"/> * Silicone only available in Beige * Brown and Black are only available in 20-30 mmHg *Brown and Black are only available in LONG length

JOBST® Bella™ Lite Armsleeve

Size	Wrist	Mid-Lower	Mid-Upper
S	15-18.9 cm (5 <sup>7/8</sup> " - 7 <sup>3/8</sup> " <sub>n</sub> )	17.8-23.2 cm (7" - 9 <sup>1/8</sup> " <sub>n</sub> )	21.5-39.5 cm (8 <sup>1/2</sup> " - 15 <sup>1/2</sup> " <sub>n</sub> )
M	19-22 cm (7 <sup>1/2</sup> " - 8 <sup>5/8</sup> " <sub>n</sub> )	23.3-28.7 cm (9 <sup>3/4</sup> " - 11 <sup>1/4</sup> " <sub>n</sub> )	27-45 cm (10 <sup>5/8</sup> " - 17 <sup>3/4</sup> " <sub>n</sub> )
L	22.1-26 cm (8 <sup>3/4</sup> " - 10 <sup>1/4</sup> " <sub>n</sub> )	28.8-34.3 cm (11 <sup>3/8</sup> " - 13 <sup>1/2</sup> " <sub>n</sub> )	32.5-50.5 cm (12 <sup>3/4</sup> " - 19 <sup>7/8</sup> " <sub>n</sub> )

Length: Regular: ≤43 cm (17") / Long: >43 cm (17")

JOBST® Bella™ Strong Armsleeve

Size	Wrist	Mid-Lower	Mid-Upper
1	14-18 cm (5 <sup>1/2</sup> " - 7")	17-21 cm (6 <sup>3/4</sup> " - 8 <sup>1/4</sup> " <sub>n</sub> )	22-31 cm (8 <sup>1/2</sup> " - 12 <sup>1/4</sup> " <sub>n</sub> )
2	14-18 cm (5 <sup>1/2</sup> " - 7")	17-21 cm (6 <sup>3/4</sup> " - 8 <sup>1/4</sup> " <sub>n</sub> )	31-41 cm (12 <sup>1/4</sup> " - 16 <sup>1/4</sup> " <sub>n</sub> )
3	14-18 cm (5 <sup>1/2</sup> " - 7")	21-25 cm (8 <sup>1/4</sup> " - 9 <sup>3/4</sup> " <sub>n</sub> )	25-34 cm (9 <sup>3/4</sup> " - 13 <sup>1/4</sup> " <sub>n</sub> )
4	16-19 cm (6 <sup>1/4</sup> " - 7 <sup>1/2</sup> " <sub>n</sub> )	21-25 cm (8 <sup>1/4</sup> " - 9 <sup>3/4</sup> " <sub>n</sub> )	28-38 cm (11" - 15")
5	17-20 cm (6 <sup>3/4</sup> " - 7 <sup>1/2</sup> " <sub>n</sub> )	21-25 cm (8 <sup>1/4</sup> " - 9 <sup>3/4</sup> " <sub>n</sub> )	34-43 cm (13 <sup>1/4</sup> " - 17")
6	14-18 cm (5 <sup>1/2</sup> " - 7")	25-29 cm (9 <sup>3/4</sup> " - 11 <sup>1/2</sup> " <sub>n</sub> )	29-38 cm (11 <sup>1/2</sup> " - 15")
7	17-20 cm (6 <sup>3/4</sup> " - 7 <sup>1/2</sup> " <sub>n</sub> )	25-29 cm (9 <sup>3/4</sup> " - 11 <sup>1/2</sup> " <sub>n</sub> )	33-43 cm (13" - 17")
8	18-21 cm (7" - 8 <sup>1/4</sup> " <sub>n</sub> )	25-29 cm (9 <sup>3/4</sup> " - 11 <sup>1/2</sup> " <sub>n</sub> )	38-48 cm (15" - 19")
9	16-19 cm (6 <sup>1/4</sup> " - 7 <sup>1/2</sup> " <sub>n</sub> )	29-33 cm (11 <sup>1/2</sup> " - 13")	32-42cm (12 <sup>1/2</sup> " - 16 <sup>1/2</sup> " <sub>n</sub> )
10	19-22 cm (7 <sup>1/2</sup> " - 8 <sup>1/2</sup> " <sub>n</sub> )	29-33 cm (11 <sup>1/2</sup> " - 13")	36-46 cm (14 <sup>1/4</sup> " - 18")

Length: Regular: ≤43 cm (17") / Long: >43 cm (17")