




Date: ___/___/___ Patient Name: _____ Therapist Name: _____

Ship to: _____



FarrowHybrid® AD Foot Compression is part compression stocking (over the foot and ankle) and part sock (ankle to knee) and may be used in place of our footpieces with our legpiece. Easy to apply and remove, they are comfortable and readily fit in patients' shoes.

FarrowHybrid®	+#SKU	Small	Medium	Large	XLarge
Male Shoe Size		3 – 6	7 – 9	10 – 12	13 – 15
Female Shoe Size		4 – 7	8 – 10	11 – 13	14 – 16
FarrowHybrid® ADI	FLFH-ADI	Qty	Qty	Qty	Qty
FarrowHybrid® ADI Wide	FLFH-ADIW	Qty	Qty	Qty	Qty
FarrowHybrid® ADIIS	FLFH-ADIIS	Qty	Qty	Qty	Qty
FarrowHybrid® ADIIS Wide	FLFH-ADIISW	Qty	Qty	Qty	Qty

Color Options	
Taupe	<input type="checkbox"/>
Black	<input type="checkbox"/>

Black available only in FarrowHybrid ADI Wide. If no color is selected, Taupe is the default option.

FarrowHybrid® Foot Compression Size Selector		
	20 – 30mmHg	30 – 40mmHg
A ^A Circumference 20 – 27cm (standard)	FLFH-ADI	FLFH-ADIIS
A ^A Circumference 28 – 35cm (wide)	FLFH-ADIW	FLFH-ADIISW

FarrowHybrids cover a wide range of circumferences in standard and wide sizes. The maximum leg circumference must be less than 60cm. FarrowHybrid® ADI and ADIW provide Class I (20-30mmHg) compression and FarrowHybrid® ADIIS and ADIISW provide Class II (30-40mmHg) compression. The FarrowHybrid® ADIIS foot compression CONTAINS SILVER and should not be used on patients allergic to silver.

FarrowHybrid® AD Foot Compression

The **FarrowHybrid® AD Foot Compression** garments are easily stretchable in the lower leg portion, which makes them easy to apply and remove, while being low profile enough to fit into most patients' current shoes. The Black color option is only available for the FarrowHybrid ADI Wide sizes.

FarrowHybrid® Measuring Instructions

Measure the circumference at the midfoot and widest calf. Midfoot circumference ranges are between 20 – 27cm for standard sizes and 28 – 35cm for the wide sizes. The widest calf circumference must be less than 60 cm and shoe size must be between 3 – 16 US. In the blank white box indicate the quantity of FarrowHybrid's desired and if applicable, the color. FarrowHybrid® AD are sold as pairs and are designed for use with Farrow brand leg garments.