



mediven® custom-made circular knit stockings

Adaptive
Together We Can Adapt to Challenges in Life

Fax (860)323-8002
Phone (860)207-8268



Date: ___/___/___

Patient Name: _____

Fitter/Therapist Name: _____

Ship to: _____

MATERIAL, COLOR, COMPRESSION						OPTIONS	
mediven® comfort color <input type="checkbox"/> wheat <input type="checkbox"/> sandstone <input type="checkbox"/> natural <input type="checkbox"/> mocha <input type="checkbox"/> navy <input type="checkbox"/> ebony style <input type="checkbox"/> calf <input type="checkbox"/> thigh <input type="checkbox"/> panty <input type="checkbox"/> maternity panty		mediven® plus color <input type="checkbox"/> beige <input type="checkbox"/> black style <input type="checkbox"/> calf <input type="checkbox"/> thigh <input type="checkbox"/> thigh w/waist <input type="checkbox"/> panty <input type="checkbox"/> maternity panty		mediven® forte color <input type="checkbox"/> caramel <input type="checkbox"/> black style <input type="checkbox"/> calf <input type="checkbox"/> thigh <input type="checkbox"/> thigh w/waist <input type="checkbox"/> panty		compressive panty* <input type="checkbox"/> 20-30 <input type="checkbox"/> 30-40 <input type="checkbox"/> 40-50 <small>*Panty compression may not be greater than legs. Compressive panty is not available in comfort.</small> silicone top band mediven® comfort: beaded (A-D) <input type="checkbox"/> 2.5cm <input type="checkbox"/> 5cm beaded (A-G) <input type="checkbox"/> 5cm lace (A-G) <input type="checkbox"/> 5cm mediven® plus & forte: beaded <input type="checkbox"/> 2.5cm <input type="checkbox"/> 5cm	
compression <input type="checkbox"/> 20-30 <input type="checkbox"/> 30-40		compression <input type="checkbox"/> 20-30 <input type="checkbox"/> 30-40 <input type="checkbox"/> 40-50		compression <input type="checkbox"/> 30-40 <input type="checkbox"/> 40-50		qty <input type="checkbox"/> left <input type="checkbox"/> right <input type="checkbox"/> pairs	
qty <input type="checkbox"/> left <input type="checkbox"/> right <input type="checkbox"/> pairs		qty <input type="checkbox"/> left <input type="checkbox"/> right <input type="checkbox"/> pairs		qty <input type="checkbox"/> left <input type="checkbox"/> right <input type="checkbox"/> pairs			
LEFT LEG circumference	LEFT LEG length to floor	WHERE TO MEASURE		RIGHT LEG circumference	RIGHT LEG length to floor	PANTY TOP length	KEY FOR CHART
							Height measurement is from each marked body location to floor LK1T Measurement from pubic bone to top of garment along the anatomical contour LK2T Measurement from base of the gluteal fold to top of garment along the anatomical contour t Measurement at waist h Measurement just above pelvic bone k Measurement at top of widest part of hips g Measurement at top of thigh at gluteal fold f Measurement at mid thigh e Measurement slightly above knee d Measurement slightly below knee c Measurement at widest part of calf b1 Measurement between ankle and waist part of calf b Measurement just above ankle bone a Measurement circumference of ball of foot y Measurement diagonally around heel over widest part of top of ankle z Measurement from heel to toe for Closed-toe stockings (enter below) or from heel to ball of foot for Open-toe stockings
7							
g	cm	g	cm	g	cm		
f	cm	f	cm	f	cm		
e	cm	e	cm	e	cm		
d	cm	d	cm	d	cm		
c	cm	c	cm	c	cm		
b1	cm	b1	cm	b1	cm		
b	cm	b	cm	b	cm		
a	cm			a	cm		
y	cm			y	cm		